



Features:

- Liquid Nitrogen Level/Temperature Alarm and Monitor
- Monitors and displays actual temperatures from -250°C to +400°C
- Type "T" thermocouple sensing device (compatible with existing sensors)
- 2-line x 20-character supertwist display
- Additional customized messaging display capability
- NO/NC Alarm Contacts
- Thermocouple Extension connectors included

Options:

- Cellular & Internet Notification System Compatible!
- Real-time clock with time/date display
- Battery backup
- Serial communication port for external data collection (non-USB)
- Monitors 1 to 8 channels
- Shorted, open or defective thermocouple sensing and alarm