

DIGITAL INSTRUMENTS

DI-120D

The Infrared Thermometer

- Quick Measurement Time (1 sec)
- Accurate Readings (+/-0.2° C)



Digital Instruments Inc.
Bulk Sales Available



FDA APPROVED



AUTO-OFF FUNCTION
Shut down automatically



FEVER ALARM
Intelligent fever warning



reddot design award

Monitor with an
award-winning
design

DI-120D



Four Color Display



Product Specification

MODEL	DI120D
BEST USE	QUICK NON-CONTACT TEMPERATURE MEASUREMENTS FOR ADULTS AND CHILDS
MEASURE RANGE	32°C~42.9°C (89.6°F~109.2°F)
ACCURACY	±0.2°C/±0.4°F
RESOLUTION	0.1°C/0.1°F
MEASURE TIME	1 SECOND FAST READING DIGITAL INFRARED THERMOMETER
WORKING TEMPERATURE	10°C~40°C; HUMIDITY: <80%
STORAGE TEMPERATURE	-25°C~+55°C; HUMIDITY: <95%
BATTERY	2*AAA , 3000 MEASUREMENT BATTERY LIFE
IDEAL MEASUREMENT DIST.	~2 INCHES (CAN BE FURTHER RISK ACCURACY)
ACCESSIBILITY	ONE MEASUREMENT BUTTON
AUDIO & VISUAL	FEVER WARNING LEDS+SOUND
UNITS	CELSIUS/FAHRENHEIT
CERTIFICATIONS/AWARD	FDA, CE, SHAPE PATENT, PRODUCT INVENTION PATENT, IF DESIGN AWARD
APPLICABLE REGULATIONS	ASTM E-1965-98, ISO80601-2-56, ISO13458, ROHS, IP22

Supplied by Digital Instruments Inc.
Custom Electronics Manufactured in
the USA for over 45 years

