

# DIGITAL INSTRUMENTS, Inc.

★ Thomas Verified

✳ COVID-19 Response

## Production and Supply Update

### DI-120D Infrared Thermometer

Digital Instruments, Inc. is proud to supply its customers bulk orders of the design award winning Infrared Thermometer.

In response to the COVID-19 crisis DI has committed to supplying this FDA approved Infrared Thermometer, manufactured with quality to control and prevent the spread of the virus.

Deemed an essential business, Digital Instruments, Inc. has remained open continuing production of all its products from its build to print assemblies to its temperature monitors for blood storage, liquid nitrogen monitors for IVF Labs and custom temperature probes.

Everyday we will continue to strive forward and support our customers and their employees.

For more information on how to secure your order of these Infrared Thermometers contact us today.

In Stock: 6,600 pc as of 6/9/2020  
Bulk Sales  
1-100k pc Pricing upon request

Contact  
sales@digitalinstruments.com  
PH: 716-874-5848 Thank You!

**-Quick Measurement Time (1 sec)**  
**-Accurate Readings (+/-0.2° C)**



  
**FDA APPROVED**

  
**AUTO-OFF FUNCTION**  
Shut down automatically

  
**FEVER ALARM**  
Intelligent fever warning



**Monitor with an award winning design**  
reddot design award



### Four Color Display



 Standby
 Normal Range
 Pay Attention
 Fever

## Bulk Sales Available

| DI-120D Tech Specs      |  |
|-------------------------|--|
| MODEL                   | DI120D   |
| BEST USE                | QUICK NON-CONTACT TEMPERATURE MEASUREMENTS FOR ADULTS AND CHILDS |
| MEASURE RANGE           | 32°C--42.9°C (89.6°F--109.2°F)                                   |
| ACCURACY                | ±0.2°C/±0.4°F  |
| RESOLUTION              | 0.1°C/0.1°F  |
| MEASURE TIME            | 1 SECOND FAST READING DIGITAL INFRARED THERMOMETER               |
| WORKING TEMPERATURE     | 10°C--40°C; HUMIDITY: <80%                                       |
| STORAGE TEMPERATURE     | -25°C--+55°C; HUMIDITY: <95%                                     |
| BATTERY                 | 2*AAA , 3000 MEASUREMENT BATTERY LIFE                            |
| IDEAL MEASUREMENT DIST. | ~2 INCHES (CAN BE FURTHER RISK ACCURACY)                         |
| ACCESSIBILITY           | ONE MEASUREMENT BUTTON   |
| AUDIO & VISUAL          | FEVER WARNING LEDS+SOUND   |
| UNITS                   | CELSIUS/FAHRENHEIT   |
| CERTIFICATIONS/AWARD    | FDA, CE, SHAPE PATENT, PRODUCT INVENTION PATENT, IF DESIGN AWARD |
| APPLICABLE REGULATIONS  | ASTM E-1965-98, ISO80601-2-56, ISO13458, ROHS, IP22              |